

# UXT antibody

Catalog Number: 116700

## Product name

UXT antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

UXT Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in HeLa cells, mouse brain tissue. Observed molecular weight by Western-blot: 18, 45 kDa

## Preparation

This antibody was obtained by immunization of UXT recombinant protein (Accession Number: NM\_004182). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP

## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116700(UXT antibody) at dilution of 1:400



IP Result of anti-UXT (IP:Catalog No:116700, 3ug; Detection:Catalog No:116700 1:300) with mouse brain tissue lysate 4000ug.