

Anti-CDH8 antibody



Catalog Number: 142010

Product name

Anti-CDH8 antibody

Immunogen

[Rat Cadherin-8 / CDH8 recombinant protein](#)

Specificity

Rat Cadherin-8 / CDH8

Antibody description

Rabbit Polyclonal to Rat CDH8

Preparation

Produced in rabbits immunized with purified, recombinant Rat Cadherin-8 / CDH8 . Cadherin-8 / CDH8 specific IgG was purified by Rat Cadherin-8 / CDH8 affinity chromatography.

Formulation

0.2 µm filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC-P


Dilutions

ELISA: 0.1-0.2 µg/mL

This antibody can be used at 0.1-0.2 µg/mL with the appropriate secondary reagents to detect Rat Cadherin-8 / CDH8.

IHC-P: 0.1-2 µg/mL

Validations



Rat Cadherin-8 / CDH8
Immunohistochemistry(IHC) 15017