

# NA14 antibody

Catalog Number: 113010

## Product name

NA14 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

NA14 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer. Positive WB detected in mouse testis tissue, mouse thymus tissue, human testis tissue. Observed molecular weight by Western-blot: 14kd

## Preparation

This antibody was obtained by immunization of NA14 recombinant protein (Accession Number: NM\_003731). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC

## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:113010(NA14 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:113010(NA14 antibody) at dilution of 1:50 (under 10x lens)