

Anti-PROS1 antibody

Catalog Number: 175010

Product name

Anti-PROS1 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to PROS1

Preparation

This antigen of this antibody was synthetic peptide within c-terminal human vitamin k-dependent protein s.

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:1,000

Validations




Fig1: Western blot analysis on human serum using anti- Vitamin K-dependent protein S polyclonal antibody.