

# ZNF600-Specific antibody



Catalog Number: 117010

## Product name

ZNF600-Specific antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

ZNF600-Specific Rabbit Polyclonal antibody.  
Positive WB detected in mouse kidney tissue.  
Observed molecular weight by Western-blot: 83 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_198457).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:117010(ZNF600-Specific Antibody) at dilution of 1:300