

# CD23 antibody

Catalog Number: 109010

## Product name

CD23 antibody

## Specificity

Human; other species not tested.

## Antibody description

CD23 Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue, human lymphoma tissue. Positive FC detected in Raji cells. Positive WB detected in Raji cells. Observed molecular weight by Western-blot: 38-42 kDa

## Preparation

This antibody was obtained by immunization of CD23 recombinant protein (Accession Number: NM\_001220500). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC, FC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:109010(FCER2,CD23 Antibody) at dilution of



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:109010(FCER2,CD23 Antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:109010(FCER2,CD23 Antibody) at dilution of

# CD23 antibody



Catalog Number: 109010

1:50 (under 40x lens)



1X10<sup>6</sup> Raji cells were stained with 0.2ug FCER2,CD23 antibody (Catalog No:109010, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.