

# RAB11FIP5 antibody



Catalog Number: 114410

## Product name

RAB11FIP5 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

RAB11FIP5 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, HEK-293 cells, human heart tissue, mouse lung tissue. Observed molecular weight by Western-blot: 66-72 kDa

## Preparation

This antibody was obtained by immunization of RAB11FIP5 recombinant protein (Accession Number: NM\_015470). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:114410(RAB11FIP5 antibody) at dilution of 1:500