

Anti-APLN antibody

Catalog Number: 176410

Product name

Anti-APLN antibody

Specificity

Human, Mouse, Rat

Antibody description

Rabbit polyclonal antibody to APLN

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from the c-terminal region of human apln. aa range:28-77

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC-P

Dilutions

WB: 1:500-1:2000

IHC-P: 1:100-300

Validations




Fig1:: Western Blot analysis of 22RV-1, H460, U87, mouse brain cells using Apelin Polyclonal Antibody.. Secondary antibody (catalog#: HA1001) was diluted at 1:20000




Fig2:: Immunohistochemical analysis of paraffin-embedded rat-brain, antibody was diluted at 1:100