

TRPM1 antibody



Catalog Number: 116410

Product name

TRPM1 antibody

Specificity

Human; other species not tested.

Antibody description

TRPM1 Rabbit Polyclonal antibody. Positive IF detected in SH-SY5Y cells.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002420).
Purification method: Antigen affinity purified.

Formulation

0.1M NaHCO₃, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IF

Dilutions

Recommended Dilution:

IF: 1:10-1:100

Validations



Immunofluorescent analysis of SH-SY5Y cells using Catalog No:116410(TRPM1 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG