

BAGE4 antibody

Catalog Number: 108410

Product name

BAGE4 antibody

Specificity

Human; other species not tested.

Antibody description

BAGE4 Rabbit Polyclonal antibody. Positive IHC detected in human lung cancer tissue, human testis tissue.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_181704).
Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:108410(BAGE4 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:108410(BAGE4 antibody) at dilution of 1:50 (under 40x lens)