

# NELFB antibody

Catalog Number: 113110

## Product name

NELFB antibody

IP: 1:200-1:2000

IF: 1:10-1:100

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

NELFB Rabbit Polyclonal antibody. Positive IF detected in BxPC-3 cells. Positive WB detected in BxPC-3 cells, HepG2 cells. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 66 kDa

## Preparation

This antibody was obtained by immunization of NELFB recombinant protein (Accession Number: BC011892). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IP, IF, WB

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Validations




BxPC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:113110(NELF-B Antibody) at dilution of



1:1000

IP Result of anti-NELF-B (IP:Catalog No:113110, 3ug; Detection:Catalog No:113110 1:500) with mouse brain tissue lysate 6000ug.



Immunofluorescent analysis of BxPC-3 cells using

# NELFB antibody



Catalog Number: 113110

---

Catalog No:113110(NELF-B Antibody) at dilution of 1:25 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)