

SENP5-Specific antibody



Catalog Number: 115110

Product name

SENP5-Specific antibody

IP: 1:200-1:1000

IHC: 1:20-1:200

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SENP5-Specific Rabbit Polyclonal antibody. Positive IHC detected in human cervical cancer tissue, human placenta tissue. Positive WB detected in HeLa cells, HepG2 cells, mouse liver tissue. Positive IP detected in HeLa cells. Observed molecular weight by Western-blot: 97 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_152699). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115110(SENP5-Specific antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human cervical cancer using Catalog No:115110(SENP5-Specific antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded

SENP5-Specific antibody

Catalog Number: 115110

human cervical cancer using Catalog No:115110(SENP5-Specific antibody) at dilution of 1:100 (under 40x lens)

IP Result of anti-SENP5-Specific (IP:Catalog No:115110, 4ug; Detection:Catalog No:115110 1:300) with HeLa cells lysate 520ug.

