

Mogamulizumab reference antibody



Catalog Number: 177110

Product name

Mogamulizumab reference antibody

Specificity

Human CCR4 / CD194

Antibody description

Anti-CCR4 / CD194 Antibody (Mogamulizumab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

1×PBS, pH 6.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG

Applications


ELISA

Validations



SDS-PAGE

Anti-CCR4 / CD194 Antibody (mogamulizumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 90%.



SEC-HPLC

The purity of Anti-CCR4 / CD194 Antibody (mogamulizumab) is more than 95%, determined by SEC-HPLC.