

Anti-Phostensin (4B12) Mouse antibody



Catalog Number: 168110

Product name

Anti-Phostensin (4B12) Mouse antibody

Specificity

Human

Antibody description

Phostensin (4B12) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human Phostensin protein fragments expressed in E.coli.

Formulation

Ascites containing with 0.03% Proclin300 and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG1


Applications

WB

Dilutions

WB: 1/200

Validations



Western blot detection of Phostensin in HeLa, Raji, Ramos and Leucocyte cell lysates and using Phostensin mouse mAb (1:200 diluted). Predicted band size: 120 kDa. Observed band size: 120 kDa.