

# ANXA2R antibody



Catalog Number: 108110

## Product name

ANXA2R antibody

## Specificity

Human; other species not tested.

## Antibody description

ANXA2R Rabbit Polyclonal antibody. Positive IHC detected in human prostate cancer tissue, human brain tissue.

## Preparation

This antibody was obtained by immunization of ANXA2R recombinant protein (Accession Number: NM\_001382352). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:20-1:200

## Validations



Immunohistochemistry of paraffin-embedded human prostate cancer tissue slide using Catalog No:108110(C5orf39 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human prostate cancer tissue slide using Catalog No:108110(C5orf39 Antibody) at dilution of 1:50 (under 40x lens)