

Mouse Alpha 1 Antitrypsin / Serpin A1a (Fc Tag) recombinant protein



Catalog Number: 504510

General Information

Protein Construction

A DNA sequence encoding the mouse Serpina1a (XP_003945711.1) (Met1-Lys413) was expressed with the Fc region of human IgG1 at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Activity

Measured by its ability to inhibit trypsin cleavage of a fluorogenic peptide substrate, Mca-Arg-Pro-Lys-Pro-Val-Glu-Nval-Trp-Arg-Lys (DNP)-NH₂. The IC₅₀ value is < 1 nM.

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Glu 25

Molecular Mass

The recombinant mouse Serpina1a comprises 627 amino acids and has a predicted molecular mass of 70.2 kDa. The apparent molecular mass of the protein is approximately 67-77 kDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Mouse Alpha 1 Antitrypsin / SerpinA1a Protein (Fc Tag) SDS-PAGE