

Human TPM4 / tropomyosin 4 (His Tag) recombinant protein



Catalog Number: 500910

General Information

Gene Name Synonym

TM30p1; Tropomyosin-4

Protein Construction

A DNA sequence encoding the human TPM4 (P67936) (Ala2-Ile248) was expressed with an N-terminal polyhistidine tag.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per µg of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant human TPM4 comprises 267

amino acids and has a predicted molecular mass of 30.8 kDa. The apparent molecular mass of the protein is approximately 35-40 kDa in SDS-PAGE under reducing conditions due to glycosylation.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human TPM4 / Tropomyosin 4 Protein (His Tag)
SDS-PAGE