

GCN5 antibody



Catalog Number: 110910

Product name

GCN5 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

GCN5 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, mouse ovary tissue. Observed molecular weight by Western-blot: 94 kDa

Preparation

This antibody was obtained by immunization of GCN5 recombinant protein (Accession Number: NM_021078). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:110910(KAT2A antibody) at dilution of 1:800