

KANK2 antibody

Catalog Number: 111910

Product name

KANK2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

KANK2 Rabbit Polyclonal antibody. Positive IF detected in Hela cells. Positive IHC detected in human colon tissue, human kidney tissue.

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001136191). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, IF


Dilutions

Recommended Dilution:


IHC: 1:20-1:200

IF: 1:10-1:100


Validations



Immunohistochemical of paraffin-embedded human colon using Catalog No:111910(KANK2 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human colon using Catalog No:111910(KANK2 antibody) at dilution of 1:100 (under 40x lens)



Immunofluorescent analysis of Hela cells, using KANK2 antibody Catalog No:111910 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).