

AFM antibody

Catalog Number: 107910

Product name

AFM antibody

Specificity

Human; other species not tested.

Antibody description

AFM Rabbit Polyclonal antibody. Positive WB detected in human blood tissue. Positive IP detected in human plasma tissue. Observed molecular weight by Western-blot: 70-85 kDa

Preparation

This antibody was obtained by immunization of AFM recombinant protein (Accession Number: NM_001133). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:1000-1:10000

Validations



human blood tissue were subjected to SDS PAGE followed by western blot with Catalog No:107910(AFMM Antibody) at dilution of 1:2000



IP Result of anti-AFM (IP:Catalog No:107910, 4ug; Detection:Catalog No:107910 1:2000) with human plasma tissue lysate 4000ug.