

Human UCHL1 (His tag) recombinant protein



Catalog Number: 519910

General Information

Gene Name Synonym

Ubiquitin Carboxyl-Terminal Hydrolase Isozyme L1, UCH-L1, Neuron Cytoplasmic Protein 9.5, PGP 9.5, PGP9.5, Ubiquitin Thioesterase L1, UCHL1

Protein Construction

A DNA sequence encoding the Human UCHL1 Met1-Ala223 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 250mM NaCl, 1mM DTT, 10% Glycerol, pH 7.5.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human UCHL1 (His tag) recombinant protein