

# EFHD2 antibody



Catalog Number: 110210

## Product name

EFHD2 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

EFHD2 Rabbit Polyclonal antibody. Positive IHC detected in human medulloblastoma tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 27 kDa

## Preparation

This antibody was obtained by immunization of EFHD2 recombinant protein (Accession Number: BC013648). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

## Validations



WB result of anti-EFHD2 (Catalog No:110210, 1:200) with mouse brain tissue lysate by Dr. Wu Zhiyuan.



Immunohistochemical of paraffin-embedded human medulloblastoma using Catalog No:110210(EFHD2 antibody) at dilution of 1:50 (under 10x lens)