

Human Tryptophan Hydroxylase 1/TPH1 (His Tag) recombinant protein



Catalog Number: 502210

General Information

Gene Name Synonym

Tryptophan 5-monoxygenase 1

Protein Construction

A DNA sequence encoding the human TPH1 (P17752-1) (Ile 2-Ile 444) was expressed, with a polyhistidine tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human TPH1 consisting of 458 amino acids and has a calculated molecular mass

of 52.7 kDa. It migrates as an approximately 48 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile 20mM Tris, 200mM NaCl, 10% glycerol, pH 8.0

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human Tryptophan Hydroxylase 1 / TPH1 Protein (His Tag) SDS-PAGE