

IKBKE Polyclonal Antibody



Catalog Number: 163210

Product name

IKBKE Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to IKBKE

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 460-620 of human IKBKE (NP_001180251.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IP


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:200

IP 1:50 - 1:200

Validations



Immunoprecipitation - IKBKE Polyclonal Antibody

Immunoprecipitation analysis of 200ug extracts of K-562 cells, using 3 ug IKBKE antibody . Western blot was performed from the immunoprecipitate using IKBKE antibody at a dilution of 1:1000.