

Prnd antibody

Catalog Number: 114210

Product name

Prnd antibody

Specificity

human, mouse; other species not tested.

Antibody description

Prnd Rabbit Polyclonal antibody. Positive WB detected in human skeletal muscle tissue. Observed molecular weight by Western-blot: 34kd

Preparation

This antibody was obtained by immunization of Prnd recombinant protein (Accession Number: NM_001278520). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations

human skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:114210(Prnd antibody) at dilution of 1:1000