

Ecromeximab reference antibody



Catalog Number: 177210

Product name

IgG1

Ecromeximab reference antibody

Applications

Specificity


ELISA

Human Disialoganglioside GD3

Validations

Antibody description

Anti-Disialoganglioside GD3 Antibody
(Ecromeximab)



Preparation


Recombinant expression and purified from CHO cells.

SDS-PAGE

Anti-GD3 Antibody (ecromeximab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0



Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

SEC-HPLC

The purity of Anti-GD3 Antibody (ecromeximab) is 96.37%, determined by SEC-HPLC.

Clonality

Monoclonal

Ig Type