

Anti-Human Serum Albumin (5G11) Mouse antibody



Catalog Number: 168210

Product name

Monoclonal

Anti-Human Serum Albumin (5G11) Mouse antibody

Ig Type

Mouse IgG1

Specificity

Human


Applications

LFIA, ELISA

Antibody description

Human Serum Albumin (5G11) Mouse monoclonal antibody

Validations



Preparation

Antigen: Native Human Serum Albumin..

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300.

Storage

Store at 4°C short term.Avoid freeze / thaw cycle.

Lateral Flow Immunoassay(LFIA) assay for anti-Human Serum Albumin mouse mAb (competitive assay).

Clonality