

# Anti-ZYX Rabbit antibody



Catalog Number: 168610

## Product name

Anti-ZYX Rabbit antibody

## Specificity

Human, Mouse

## Antibody description

ZYX Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human Zyxin

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, ICC/IF, FC, IP

## Dilutions

WB: 1/2000-1/10000


IHC: 1/50

ICC/IF: 1/20-1/100

FC: 1/200-1/1000

IP: 1/20

## Validations



Western blot detection of ZYX in C6,Rat Brain lysates using ZYX Rabbit mAb(1:500 diluted).Predicted band size:61kDa.Observed band size:78kDa.