

Human COASY / CoA synthase (His Tag) recombinant protein



Catalog Number: 501310

General Information

Gene Name Synonym

NBP; POV-2; Phosphopantetheine adenylyltransferase; Dephospho-CoA pyrophosphorylase; Pantetheine-phosphate adenylyltransferase; Dephospho-CoA kinase; Dephosphocoenzyme A kinase

Protein Construction

A DNA sequence encoding the human COASY (Q13057-1) (Met296-Asp564) was expressed with a polyhistidine tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

> 90 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

His

Molecular Mass

The recombinant human COASY consists of 287 amino acids and predicts a molecular mass of 32.5 KDa. It migrates as an approximately 33 KDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human COASY / CoA synthase Protein (His Tag)
SDS-PAGE