

NOV antibody



Catalog Number: 113310

Product name

NOV antibody

Recommended Dilution:

WB: 1:200-1:2000

Specificity

Human, Mouse, Rat; other species not tested.

IHC: 1:20-1:200


Validations



Antibody description

NOV Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue. Positive IHC detected in human nephroblastoma tissue, human prostate cancer tissue. Observed molecular weight by Western-blot: 45kd


mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:113310(NOV antibody) at dilution of 1:500



Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_002514). Purification method: Antigen affinity purified.

Immunohistochemical of paraffin-embedded human nephroblastoma using Catalog No:113310(NOV antibody) at dilution of 1:50 (under 10x lens)



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Immunohistochemical of paraffin-embedded human nephroblastoma using Catalog No:113310(NOV antibody) at dilution of 1:50 (under 40x lens)

Dilutions