

Phospho-AKT2-S474 pAb



Catalog Number: 166310

Product name

Phospho-AKT2-S474 pAb

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to Phospho-AKT2-S474

Preparation

Antigen: A synthetic phosphorylated peptide around S474 of human AKT2 (NP_001617.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - Phospho-AKT2-S474 pAb

Western blot analysis of extracts of 293 cells, using Phospho-AKT2-S474 antibody 1:1000 dilution..Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% BSA.