

Human POMGNT1/MGAT1.2/UNQ746/PRO1475 (His tag) recombinant protein



Catalog Number: 519310

General Information

Gene Name Synonym

Protein O-Linked-Mannose Beta-1,2-N-Acetylglucosaminyltransferase 1, POMGnT1, UDP-GlcNAc:Alpha-D-Mannoside Beta-1,2-N-Acetylglucosaminyltransferase I.2, GnT I.2, POMGNT1, MGAT1.2

Protein Construction

A DNA sequence encoding the Human POMGNT1/MGAT1.2/UNQ746/PRO1475 Leu59-Thr660 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM TrisHCl, 150mM NaCl, 10% Glycerol, pH 8.5.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human POMGNT1/MGAT1.2/UNQ746/PRO1475 (His tag) recombinant protein