

# FOXP3 antibody



Catalog Number: 110710

---

## Product name

FOXP3 antibody

## Immunogen

[Human FOXP3 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse; other species not tested.

## Antibody description

FOXP3 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, A375 cells, HeLa cells, RAW264.7 cells. Observed molecular weight by Western-blot: 47-50 kDa

## Preparation

This antibody was obtained by immunization of FOXP3 recombinant protein (Accession Number: NM\_014009). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:110710(FOXP3 antibody) at dilution of 1:1000