

# MMS19 antibody



Catalog Number: 112710

---

## Product name

MMS19 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

MMS19 Rabbit Polyclonal antibody. Positive WB detected in human brain tissue, HeLa cells, L02 cells, mouse brain tissue, mouse kidney tissue, mouse testis tissue. Observed molecular weight by Western-blot: 103 kDa

## Preparation

This antibody was obtained by immunization of MMS19 recombinant protein (Accession Number: BC006575). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:112710(MMS19 antibody) at dilution of 1:500