

Human GITR / TNFRSF18 (His Tag) recombinant protein



Catalog Number: 503710

General Information

Gene Name Synonym

Activation-inducible TNFR family receptor;
Glucocorticoid-induced TNFR-related protein

Protein Construction

A DNA sequence encoding the human TNFRSF18 (NP_683700.1) (Met1-Glu161) was expressed with a polyhistidine tag at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Purity

> 90 % as determined by SDS-PAGE.

Endotoxin

< 1.0 EU per μg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 26

Molecular Mass

The recombinant human TNFRSF18 consists 147 amino acids and predicts a molecular mass of 16 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.
1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human GITR / TNFRSF18 Protein (His Tag)