

peptide YY antibody



Catalog Number: 113710

Product name

peptide YY antibody

Specificity

Human; other species not tested.

Antibody description

peptide YY Rabbit Polyclonal antibody. Positive IHC detected in human pancreas tissue, human colon tissue, human small intestine tissue.

Preparation

This antibody was obtained by immunization of peptide YY recombinant protein (Accession Number: XM_011525035). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemistry of paraffin-embedded human pancreas slide using Catalog No:113710(PYY Antibody) at dilution of 1:50



Immunohistochemistry of paraffin-embedded human pancreas slide using Catalog No:113710(PYY Antibody) at dilution of 1:50