

Anti-ING5 antibody



Catalog Number: 104710

Product name

Anti-ING5 antibody

Immunogen

[Human ING5 \(GST Tag\) recombinant protein](#)

Specificity

Human ING5

Antibody description

Rabbit monoclonal to ING5

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human ING5 (rh ING5; Q8WYH8-2; Glu36-Gln84).

Formulation

0.2 μ m filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP

Dilutions


WB: 2-10 μ g/ml

ELISA: 0.1-0.2 μ g/ml


This antibody can be used at 0.1-0.2 μ g/ml with the appropriate secondary reagents to detect Human ING5.

IP: 1-4 μ L/mg of lysate

Validations



ING5 Antibody, Rabbit MAb, Western blot



ING5 Antibody, Rabbit MAb, Immunoprecipitation