

VAMP8 antibody

Catalog Number: 116710

Product name

VAMP8 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

VAMP8 Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells. Observed molecular weight by Western-blot: 15 kDa

Preparation

This antibody was obtained by immunization of VAMP8 recombinant protein (Accession Number: NM_003761). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:116710(VAMP8 Antibody) at dilution of 1:300