

Anti-Activin Receptor Type IIB/ACVR2B Rabbit antibody



Catalog Number: 167710

Product name

Monoclonal

Anti-Activin Receptor Type IIB/ACVR2B Rabbit antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse, Rat

Applications

WB

Antibody description

Activin Receptor Type IIB/ACVR2B Rabbit monoclonal antibody


Dilutions

WB: 1/1000

Preparation

Antigen: A synthetic peptide of human Activin Receptor Type IIB/ACVR2B

Validations



Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Western blot detection of Activin Receptor Type IIB/ACVR2B in CHO-K1,mouse Brain using Activin Receptor Type IIB/ACVR2B Rabbit mAb(1:1000 diluted)

Clonality