

# CYP4Z1-Specific antibody



Catalog Number: 109710

## Product name

CYP4Z1-Specific antibody

## Specificity

Human; other species not tested.

## Antibody description

CYP4Z1-Specific Rabbit Polyclonal antibody.  
Positive WB detected in MDA-MB-453s cells, HeLa cells. Observed molecular weight by Western-blot: 59kd

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_178134).  
Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol  
pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



MDA-MB-453s cells were subjected to SDS PAGE followed by western blot with Catalog No:109710(CYP4Z1-Specific antibody) at dilution of 1:300