

Anti-DCHS1 antibody



Catalog Number: 175020

Product name

Anti-DCHS1 antibody

Ig Type

Rabbit IgG

Specificity

Human

Applications

WB

Antibody description

Rabbit polyclonal antibody to DCHS1

Dilutions

WB: 1:500-1:1,000

Preparation

ICC: 1:50-1:100

This antigen of this antibody was synthetic peptide within human protocadherin-16 aa 629-673.

Formulation

Liquid, 1*PBS (pH7.4), 0.2% BSA, 40% Glycerol.
Preservative: 0.05% Sodium Azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Validations




Fig1: Western blot analysis on the recombinant protein of protocadherin-16 using anti-protocadherin-16 polyclonal antibody.