

Anti-LOC100653061 antibody



Catalog Number: 176020

Product name

Anti-LOC100653061 antibody

Specificity

Human

Antibody description

Rabbit polyclonal antibody to LOC100653061

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human gusb11. aa range:51-100

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1.; Western Blot analysis of K562 cells using GUSBP1 Polyclonal Antibody