

APLP1 antibody

Catalog Number: 108020

Product name

APLP1 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

Specificity

Human, Mouse; other species not tested.

Antibody description

APLP1 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse brain tissue. Positive IHC detected in human brain tissue. Observed molecular weight by Western-blot: 72 kDa

Preparation

This antibody was obtained by immunization of APLP1 recombinant protein (Accession Number: NM_005166). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IHC, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108020(APLP1 Antibody) at dilution of 1:1000



IP Result of anti-APLP1 (IP:Catalog No:108020, 4ug; Detection:Catalog No:108020 1:600) with mouse brain tissue lysate 3440ug.




Immunohistochemistry of paraffin-embedded human brain tissue slide using Catalog No:108020 antibody.

APLP1 antibody

Catalog Number: 108020

No:108020(APLP1 Antibody) at dilution of 1:50
(under 10x lens)



Immunohistochemistry of paraffin-embedded
human brain tissue slide using Catalog
No:108020(APLP1 Antibody) at dilution of 1:50
(under 40x lens)