

Cynomolgus TNFRSF21/DR6 (His Tag) recombinant protein



Catalog Number: 504420

General Information

Protein Construction

A DNA sequence encoding the cynomolgus TNFRSF21 (XP_005552846.1) (Met1-Leu350) was expressed with a polyhistidine tag at the C-terminus.

Organism

Cynomolgus

Expression Host

Human Cells

QC Testing

Activity

Measured by its ability to bind with APP-Fc(Cat:501404) in a functional ELISA.

Purity

> 95 % as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL method.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Gln 42

Molecular Mass

The recombinant cynomolgus TNFRSF21 consists 320 amino acids and predicts a molecular mass of 35 kDa.

Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Cynomolgus TNFRSF21 / DR6 Protein (His Tag)
SDS-PAGE