

# Mcp1 antibody



Catalog Number: 107420

## Product name

Mcp1 antibody

## Specificity

Mouse; other species not tested.

## Antibody description

Mcp1 Mouse Monoclonal antibody. Positive WB detected in LPS treated RAW 264.7 cells. Observed molecular weight by Western-blot: 16 kDa

## Preparation

This antibody was obtained by immunization of Mcp1 recombinant protein (Accession Number: NM-011333). Purification method: Protein G purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



LPS treated RAW 264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:107420(Mcp1 Antibody) at dilution of 1:1000