

Drozitumab (Synonyms: USP MAB 003, MONOCLONAL IGG1 [USP-RS]; ANTI-DR5 RHUMAB-DR5)



Catalog Number: 176820

Product name

Drozitumab (Synonyms: USP MAB 003, MONOCLONAL IGG1 [USP-RS]; ANTI-DR5 RHUMAB-DR5)

Specificity

Human TNFRSF10B / TRAILR2 / CD262

Antibody description

Anti-TNFRSF10B / TRAILR2 / CD262 Antibody (Drozitumab)

Preparation

Recombinant expression and purified from CHO cells.

Formulation

0.1 M Pro-Ac, 20 mM Arg, pH 5.0

Storage

-80°C for 2 years under sterile conditions; -20°C for 1 year under sterile conditions; Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal


Ig Type

IgG1

Applications


ELISA

Validations




SDS-PAGE

Anti-TNFRSF10B / TRAILR2 / CD262 Antibody (drozitumab) on SDS-PAGE under reducing (R) condition. The purity of the protein is greater than 95%.



SEC-HPLC

The purity of Anti-TNFRSF10B / TRAILR2 / CD262 Antibody (drozitumab) is 99.68%, determined by



SEC-HPLC.

Bioactivity: ELISA

Immobilized Cynomolgus TRAIL R2/DR5/ TNFRSF10B Protein, His Tag at 2 ug/mL can bind drozitumab, EC50=0.1257 ug/mL.