

WEE2 antibody



Catalog Number: 116820

Product name

WEE2 antibody

Specificity

Human; other species not tested.

Antibody description

WEE2 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, human ovary tissue. Observed molecular weight by Western-blot: 63-66 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001105558). Purification method: Antigen affinity purified.

Formulation

0.1M NaHCO₃, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:116820(WEE2 antibody) at dilution of 1:500