

Calcyphosine 2 antibody



Catalog Number: 108820

Product name

Calcyphosine 2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

Calcyphosine 2 Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor tissue. Positive WB detected in mouse testis tissue. Observed molecular weight by Western-blot: 54 kDa

Preparation

This antibody was obtained by immunization of Calcyphosine 2 recombinant protein (Accession Number: XM_005269194). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:108820(CAPS2 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:108820(CAPS2 antibody) at dilution of 1:50 (under 10x lens)