

# RPLP2 Polyclonal Antibody



Catalog Number: 165120

---

## Product name

RPLP2 Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to RPLP2

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-115 of human RPLP2 (NP\_000995.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IF

## Dilutions

IF 1:50 - 1:100

## Validations



Immunofluorescence - RPLP2 Polyclonal Antibody

Immunofluorescence analysis of U2OS cells using RPLP2 antibody . Blue: DAPI for nuclear staining.