

TICAM2 antibody

Catalog Number: 116120

Product name

TICAM2 antibody

Specificity

Human; other species not tested.

Antibody description

TICAM2 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, HeLa cells, human testis tissue. Positive FC detected in HeLa cells. Positive IHC detected in human testis tissue, human prostate hyperplasia tissue. Observed molecular weight by Western-blot: 32 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_021649). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, FC


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:116120(TICAM2 antibody) at dilution of 1:500



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:116120(TICAM2 Antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:116120(TICAM2 Antibody) at dilution of 1:50 (under 10x lens)

TICAM2 antibody



Catalog Number: 116120

No:116120(TICAM2 Antibody) at dilution of 1:50
(under 40x lens)



1X10⁶ HeLa cells were stained with 0.2ug TICAM2 antibody (Catalog No:116120, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.